
cub[®]

Comfortable upright birth

Information Guide

Better for mama, better for baby

cub[®] contents

About the cub[®]

Cleaning your cub[®]

Why upright positions matter

Instructions for inflating your cub[®]

Safety guidance for using the cub[®]

Physiotherapy guidance for cub[®] use if you have pelvic girdle pain in pregnancy

Positions you can use your cub[®]

Hospital or Birth Centre cleaning guide

Upright positions to use with your cub[®] during labour and birth

Troubleshooting

Research on upright birth

Information in different languages

About cub®

The cub provides additional support and comfort during late pregnancy, labour and birth.

What is the CUB?

The **cub**® is a modern innovation of the traditional birth stool and the widely used birth ball. Its innovative design supports all of the advantages of an active labour and an upright position for birth. Lightweight, adaptable, comfortable and easy to clean and can be used in any environment.

The **cub**® is a simple, easy to use invention that supports a range of upright positions, recommended by the World Health Organisation.

Upright positions can help to decrease complications that can lead to emergency interventions and even caesarean sections.

What does it do?

The **cub**® is designed to provide all of the advantages of an active labour while supporting a variety of positions for giving birth. The **cub**®'s patented, innovative design has removed the disadvantages associated with the traditional type birth stool or the common birth ball.

The **cub**® supports the needs and choices of mothers by facilitating mobility, comfort and relaxation during pregnancy, labour and birth.

Even when things are not as straight forward as planned the **cub**® can continue to be a useful, adaptable support giving the mother the best opportunity for a safe, healthy and natural birth.

Kneeling, squatting, all fours as well as sitting are all easily and comfortably achieved with the **cub**®.

Where can the **cub**® be used?

The **cub**® can be used at home, in hospitals or in a birth centre.

Lightweight and portable the **cub**® can be moved to wherever it is needed quickly and easily, helping to maintain an upright position for birth wherever you are.

Caring for your **cub**®

The **cub**® is a semi disposable product, however if cared for properly, it will be suitable for multiple uses.

Proper care includes:

Not allowing anyone to “play” with The **cub®.** This includes partners, carers, children, pets and visitors.

Do not allow more than one person to sit on the **cub® at any one time.**

The **cub® requires regular cleaning & inspection to remain in good condition.**

When not in use, the cub should be stored deflated.

Keep away from sources of heat, uneven surfaces and sharp objects.

You can find details on how to clean your **cub® further on in this brochure.**

Cleaning your cub®

The cub® should be cleaned thoroughly between uses.

We recommend chemical cleaning using one of the following cleaners:

Antibacterial spray or wipes, sporicidal spray or wipes or a diluted hypochlorite (bleach) solution.

Personal Protective Equipment (PPE), such as gloves and aprons, should be worn (where appropriate).

To clean whilst the cub® is inflated:

Rinse off all visible surface contamination with clean water until it looks clean on every surface, pay attention to seam, valve and base areas.

Dry the cub® with a clean disposable cloth, such as a household towel.

Spray or wipe the cub® with a selected cleanser from the above recommended products.

Allow the cub® to air dry completely. Deflate and store the cub®. Do not store in very warm or very cold environments as this can degrade the material.

Why are upright positions important during birth?

Research has shown that adopting upright positions during labour and birth can have many advantages for both mother and baby.

The chance of your baby becoming distressed during labour when you are upright is decreased by

54%

Being upright during labour and birth can increase the available space within the pelvis by up to

30%

Giving more room for your baby to be born.

Your chance of needing a cut in your perineum to help baby be born is reduced by

21%


If you are upright.

Your chance of needing medical assistance e.g. forceps is reduced by

23%

If you are upright during birth

Emergency caesarean sections
are reduced by

29%

If you birth in an upright position

Some studies have shown that an upright position in labour is as effective as an Oxytocin (Pitocin) drip in stimulating contractions.

The length of the early and middle parts of your labour can be shorter when adopting an upright position

The time taken to push your baby out is shorter with upright positions.

You will have less chance of needing stitches for a tear and less swelling in an upright position during birth,

Instructions for your cub[®] inflation

Make sure your cub[®] is at room temperature before inflation.

Unpack the cub[®] and unfold.

Locate the two safety inflation valves at the back of the cub[®].

Insert the nozzle of the inflation pump into the black section of the valve.

Inflate each chamber of the cub[®] until it is firm to touch.

Simply remove the black section by pulling on the clear tab on the valves and allow the air to expel.

The cub[®]'s double section valves are a SAFETY feature of the product. Once fully inflated and during use, the valves must be pushed in firmly. Your cub[®] is now ready to use.

As your cub[®] is used the material stretches slightly so you may need to top up the air.

If you are using the cub[®] Cover: Half inflate the cub[®], put the cover over, once in place finish inflating your cub[®], adjusting cover as required.

Using your cub® safely

The wide, flat base provides enhanced stability for users. The maximum user weight is 120kg.

The **cub**® should never be used to sit or lean on by more than one person at any one time.

Never stand on or try to balance on the **cub**®. Feet or knees should always be in contact with a stable surface such as the floor.

If the mother has been in water or is wet for any reason place an absorbent towel on the **cub**® before sitting/leaning over or against it.

The **cub**® should never be used as child's toy, child's seat or in a swimming pool.

Keep away from open sources of heat – the **cub**® is PVC, which can distort & melt when close to direct heat.

Keep away from potential source of puncture, particularly hard uneven surfaces or sharp objects.

Do not inflate the **cub**® unless it has been stored at room temperature for at least 4 hours.


Positions you can use your cub® in

The cub® can be adjusted to be used in a single or dual inflation as required.

Asymmetrical
(one knee)
kneeling position
(Single inflation)

Sitting upright
(Single inflation)

Forward kneeling
position on a bed
(Single inflation)

Side lying position
(Single inflation)

Supported upright
kneeling
(Dual inflation)

Sitting upright
(Dual inflation)

Supported forward kneeling
position
(Dual inflation)

Using the cub® to
support squatting
(Dual inflation)

Physiotherapy Guidance

Physiotherapy guidance if you have Pelvic Girdle Pain (PGP) in pregnancy or during childbirth. Using the cub® to ease discomfort.

Dual chamber inflation positions

Supported kneeling position

Coccyx pain
Prolapsed intervertebral disc
Cub support height: 41cm
Hip Flexion: 90° – 100°
Pain free abduction required

Sitting upright with support


Posterior pelvic pain
Prolapsed intervertebral disc
Spinal Stenosis
Cub support height: 41 cm
Hip Flexion: 90° – 100°
Abduction: 20°

Supported forward kneeling position

Posterior pelvic pain
Coccyx pain
Spondylolisthesis
Cub support height: 41 cm
Hip Flexion: 90°
Abduction: 20°

Single chamber inflation position

Side laying position

Pubic symphysis dysfunction or with Epidural
Cub support height when used in single chamber inflation: 21cm
Hip Flexion: 80°
Abduction: 20°

Lying on the comfortable side with leg supported is a more suitable and comfortable position.
J.S. Boissonault 2001

Hospital or Birth Centre cleaning guide

Suitable cleaning products to use with the **cub**[®]:

Antibacterial spray or wipes, Haz Tabs, Chlor-Clean, Actichlor Plus, Sporidical spray or wipes* Hypochlorite.

*This includes Milton[®] sterilisation fluid. Milton Fluid is made of an aqueous solution of sodium hypochlorite and 16.5% sodium chloride. The Milton Fluid that is available to buy is 2% sodium hypochlorite. PPM concentrations are as follows:

For blood contamination:

1 part Milton 2% to 2 parts water provides a 10,000ppm cleaning solution

For body fluid contamination that is not blood stained:

1 part Milton 2% to 20 parts water provides a 1000ppm cleaning solution.

The cleaning process for the **cub**[®] We recommend that a clean, disposable incontinence sheet is placed onto the **cub**[®] if the mother is sitting on it and has any PV loss that may come into direct contact with the **cub**[®]. The **cub**[®] should be thoroughly cleaned between users using the process below. Personal protective equipment (PPE) such as gloves and aprons should be worn during the cleaning process.

While the **cub**[®] is still inflated:

Rinse – Rinse off all visible surface contamination with a clean, disposable cloth using water/detergent; paying attention to the area around the seams and valves and including the base. There should be no visible contamination left.

Dry – Dry the **cub**[®] with a clean disposable cloth, such as household kitchen towel or paper hand towels.

Clean – Clean thoroughly with your chosen, locally approved cleaning product.

Dry - Allow the **cub**[®] to air dry completely before deflating and storing.

Troubleshooting for your cub®

Top up the cub® with air – the cub® is made of PVC, which stretches with repeated use, particularly in warm environments. If the cub is frequently used, we recommend you replace it every three months

If either chamber of the **cub®** is deflating please add more air and ensure that both valves are firmly pushed in. Remember it is expected that you will need to occasionally add more air if the **cub®** is left inflated for more than 24 hours.

If a puncture is suspected, re-inflate the **cub®** and locate the puncture by rubbing a small amount of washing liquid/bubble bath over the seams & valves/area of suspected damage. Push down gently on the **cub®** & identify any small bubbles appearing. Then locate the area to apply the repair patch.

A torn seam cannot be repaired and is an indication that excessive pressure or weight has been placed on the **cub®**.

As an inflatable product, the **cub®** is not designed to last indefinitely. Frequent use will naturally degrade the product sooner than occasional/ infrequent/ single user use.

Research

On upright positions for birth

Being upright during labour and birth can increase the available space within the pelvis by 28-30% giving more room for rotation and descent [6]

Uterine contractions are more effective. Some preliminary studies have shown that an upright position in labour is as effective as an Oxytocin (Pitocin) infusion in stimulating contractions [8]

There is a 54% decreased incidence of foetal heart rate abnormalities when the mother is upright [6]

The length of the first stage of labour, between 3-10 centimetres can be significantly shortened when the mother remains spontaneously upright, active and mobile [1, 7]

Upright positions reduce the duration of the second stage of labour [3, 6]

Being upright has been shown to result in a 23% reduction in assisted deliveries - more space, less distress and a shorter labour results in fewer necessary interventions [6]

A 21% reduction in rates of episiotomy in women who give birth in an upright position [4, 6]

A 29% reduction in emergency caesarean sections [7]

Improved perineal integrity, less vulvular oedema [1, 4]

Reduced need for epidural [7]

1 Terry RR, Westcott J, O'Shea L, Kelly F. (2006).

Postpartum Outcomes in Supine Delivery by Physicians vs Non-supine Delivery by Midwives. *J Am Osteopath Association*. vol. 106 no. 4, 199-202.

2 Keen R, DiFranco J, Amis D, Albers L. (2004)

Non-Supine (e.g., Upright or Side-Lying) Positions for Birth. *Journal of Perinatal Education*. 2004 spring; 13(2): 30-34.

3 Lavender T and Mlay R. (2006)

Position in the second stage of labour for women without epidural anaesthesia: RHL commentary. The WHO Reproductive Health Library; Geneva: World Health Organization. last revised: 15 December 2006.

4 Nasir A, Noorani K.J., Korejo R. (2007)

Child birth in squatting position. Department of Obstetrics and Gynaecology, Jinnah Postgraduate Medical Centre, Karachi.

5 Shermer RH. Raines DA. (1997)

Positioning during the second stage of labor: moving back to basics *JOGNN: Journal of Obstetric, Gynaecologic, and Neonatal Nursing*. Nov-Dec; 26(6): 727-34.

6 Gupta JK, Hofmeyr GJ, Smyth R. (2012).

Position in the second stage of labour for women without epidural anaesthesia. *Cochrane Database Systematic Review*. May 16;5 :CD002006 doi:10.1002/14651858, CD002006.pub3.

7 Lawrence A, Lewis L, Hofmeyr GJ, Dowswell T, Styles C. (2013)

Maternal positions and mobility during first stage labour. *Cochrane Database of Systematic Reviews Issue 2*. Art. No.:CD003934; DOI: 10.1002/14651858.CD003934.pub2.

8 Read JA, Miller FC, Paul RH. (1981).

Randomized trial of ambulation versus oxytocin for labour enhancement: a preliminary report. *Am J Obstet Gynecol*. Mar 15;139(6)669-72.

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Tel: +44 (0) 1563 851404

Email: sales@cubmama.com

f [facebook.com/CUBsupport](https://www.facebook.com/CUBsupport)

t [@CUBcass](https://twitter.com/CUBcass)

i [@CUBsupport](https://www.instagram.com/CUBsupport)

p [@CUBsupport](https://www.pinterest.com/CUBsupport)

PVC can be recycled to make other products.

Medical grade PVC with a biocompatible coating.

We recommend the **cub**[®] is replaced every 3 months if used by multiple users in a health care setting.

The **cub**[®] can be adjusted in height for safe use on a hospital bed if required.

As the **cub**[®] warms up once inflated it becomes more flexible and stretchy, so the air may need topped up if it is left inflated.

The wide, flat base provides enhanced stability for users.

Birthsparks

Birthsparks Ltd
Unit A1A, Olympic Business Park
Drybridge Road, Dundonald
KA2 9BE Scotland

Tel: +44 (0)1563 851404
Email: sales@cubmama.com
Web U.K.: www.cubmama.com
Web U.S.A.: www.the-cub-usa-com

Product specifications:
Maximum user weight: 120 kg
Product weight: 2 kg
Total weight: 4 kg

Material 0.500mm Biocompatible PVC
Double layered thickness seat and floor covering
Two separate inflation valves/Chambers
Complies with EN71 -1, 2, 3